




Is-Like Presentation

Ortho Team

Bianca, Braxton, Jonathan, Li

Discovery



Clinical Observations



No. Procedures: 9
Total Hours: **46**



No. Procedures: 11
Total Hours: **51**



No. Procedures: 9
Total Hours: **45.5**



No. Procedures: 11
Total Hours: **57**

Total Hours

200

(199.5)

THANKS TO



Dr. Andrew Stein



Dr. Li



Dr. Olszewski



Dr. Freccero



Dr. C Tannoury




Dr. Kain




Dr. T Tannoury

And all the **nurses** and **residents** – especially **Dr. Francesca Simeone**, who were all extremely helpful and indispensable to running operations


Observed Procedures



Arthroscopic Rotator Cuff Repairs




Total and Partial Shoulder, Hip, and Knee Arthroplasties (Joint Replacements)




Hand and Wrist Procedures



Internal Fixation and Repairs of Lower Limbs




Spinal Fusions



Spinal Laminectomies (Decompressions)


Problem Selection

Rotator cuff retears



Biological
Problem

Joint replacement fitting




Not a real
problem and
tight
tolerances



Spinal fusion complications




Too complex



Skeletal Traction



Skeletal Traction






Need Statement


Orthopedic surgeons need a way to apply skeletal traction to adult patients with femoral shaft fractures that maintains correct femur length throughout treatment despite changes in patient position to prevent femur displacement that can cause further damage.



Stakeholders



Design




Current Technology Limitations



Expensive!

Either portable OR
sufficient traction

No indication of how
much force is applied



Proposed Solution

1. No weights
2. Self-contained
3. Adjustable Countertraction



Custom
Traction Winch

Cable

Force Gauge

Countertraction Angle
Adjustment Brackets

Traction
Frame

Leg Cradle
and Straps



Value Proposition

Our portable skeletal traction device allows orthopedic surgeons to apply and monitor traction, maintaining correct femur length even **during patient repositioning**.

Unlike current methods, the device is **compatible with most standard hospital beds**, and continued **traction throughout scans** enables more accurate treatment plans for better outcomes.

Development





Device Requirements

	User Need		User Need
1.	Device shall apply traction to femur fracture per current treatment requirements without excessive effort	6.	The pin must align with current traction pin use and safety standards
2.	The device shall maintain the set length and position of femur throughout treatment	7.	Device shall allow for timely evaluation, monitoring, and cleaning of pin site
3.	Device must measure and relay force readings to user	8.	Device pin shall be biocompatible
4.	Device shall be radiolucent for X-rays	9.	Device pin shall be sterilizable
5.	The device must be compatible with use of most hospital beds to allow for efficient patient transport	10.	Device packaging and labeling shall be consistent with FDA standards
		11.	Device shall be composed of reusable components



User Need	Design Input
Device shall apply traction to femur fracture per current treatment requirements without excessive effort	Device must be able to apply 5-40 lbs (10-20% of patient's body weight) of traction axially
The device shall maintain the set length and position of femur throughout treatment	Force applied by the device does not deviate more than 2.5 lbs of set force applied
Device must measure and relay force readings to user	Device shall measure applied force with accuracy of ± 2 lbs
	Force readings are displayed to user
Device shall be radiolucent for X-rays	Device must be made of components that do not obstruct the fracture during X-rays
The device must be compatible with use of most hospital beds to allow for efficient patient transport	Device fits within bounds of 36 x 80 x 24 inches

Device Inputs and Outputs

Traction Frame -
6061-T6 Aluminum




Custom Traction Winch -
Ratchet and pawl
w/ worm drive

Cable -
Vinyl-coated
steel


Force Gauge -
Indicates force
to user

Digital Solidworks Prototyping Iterations


Iteration 1




Iteration 2



Iteration 3




Iteration 4



FEA of Validation of Final Model

Verification of mechanical properties


1. Aluminum Yield Tolerance
 - a. 40,000 - 70,000 PSI
2. FEA test → stress distribution
3. SN - Curve comparison (Aluminum)



At 4.8 ksi, the design is below fatigue limits (device > 10⁸ cycles)

Stress Test Analysis

Von Mises Stress (PSI)



Max Stress of 4847 PSI

Mechanical analysis → device lifespan 5+ years minimum


Physical Prototyping

Mockup


Initial Works Like

Advanced Works Like


Is Like



Custom Winch Design




Custom Winch Design



Worm Drives!

Input and
output shaft
orientation

High torque
reduction ratios
(used 20:1)




Smooth,
incremental
adjustments

Self-locking (for
lead angles $< 5^\circ$)


Compact

Risk Analysis and Mitigation


Double Locking with
Ratchet and Pawl



Removable handle to prevent
accidental/untrained adjustment



Risk Analysis and Mitigation



Include leg cradle, padding, and straps




Verification Testing - Maintaining Traction

- 19.4 lb traction force applied at 45 degree setting
- 20 lb leg stand-in attached to force gauge
- Traction force measurements taken in intervals
- Target: ≤ 2.5 lb traction decay over 24 hrs




Verification Testing - Maintaining Traction



Verification Testing - Winch Handle

- Measuring force required to turn winch handle using compression force gauge
- Target: < 5 lbs maximum force requirement at 40 lbs of traction



Future Verification - Radiolucency

- Obtain bone equivalent material with small detailed targets
- Image baseline without traction, then material in traction
- Target: Maintain $\geq 80\%$ original SNR, ≥ 3 CNR





Future Verification - Sterility

Frame

- Wipe down test on frame using standard disinfectants
- Target: $< 6.4 \mu\text{g}/\text{cm}^2$ protein residue after cleaning

Pin

- Force testing of pouch opening and cap removal, atmospheric, mechanical, vibrational testing on pouch
- Target: $2 < x < 5$ lb seal on pouch and cap, maintain SAL 10^{-6}

Validation Plan - Traction Workflow

- Assess time needed for staff to apply or take off traction to patient stand-in compared to SOC
- Target: User can set traction at 10, 20, 30, 40 lbs with ± 2 lb accuracy and turn winch with < 5 lbs of force



Validation Plan - Frame Portability

- Test whether clinical staff can transfer patient stand-in traction from bed to scanning table
- Target: Staff can complete necessary tasks pertaining to patient transport in <5 min additional time



Validation Plan - X-Ray Compatibility

- Test frame compatibility with standard femur fracture scanning workflow, measuring time needed to complete tasks
- Target: Staff can place pin through bone proxy and can visualize fracture on scans



Deployment



Femur Fractures By The Numbers

US Population in 2021
330 Million



US incidence of femoral
shaft fractures
9.5-19 per 100,000

Femoral Shaft Fractures
Treated with Traction
~23%*



*** varies widely by location**

STANDARD
WORKFLOW
\$1,744

\$420

\$54

\$70

\$900

\$300

ARBUS MEDICAL

COST TABLE

Tray Sterilization^{1,5}

Power System • Battery Pack • Traction Pin • Kirschner Bow

Capital Depreciation & Maintenance^{2,3,4,5}

Single Use Materials⁵

Prep Materials • Dressing Supplies

Staff Time Gathering & Opening^{5,6,7,8,9}

Staff Time Applying Traction⁵

WORKFLOW
WITH TrakPak®
\$1,129

\$615
SAVINGS
PER CASE

\$829
Cost of
TrakPak®



\$100

\$200



Market Size

Total Addressable Market (TAM)

10 / 100,000 cases x 330M = 33,000

Revenue/sale	\$80+\$1800
--------------	--------------------

Sales/year	33,000+330
------------	-------------------

Total	\$3.23M
-------	----------------

Serviceable Available Market (SAM)

33,000 cases x 40% = 13,200


Revenue/sale	\$80 +\$1800
--------------	---------------------

Sales/year	13,200+132
------------	-------------------

Total	\$1.29M
-------	----------------

Existing IP Space

Patent: US6467713B1
Drum-based winch cable traction
Status: Expired



Patent: US6467713B1
Cervical traction with force measurement device
Status: Active





Fig.1

Patent: US3364925A
Cables & pulleys → constant traction
Status: Expired



U.S. Patent Aug. 21, 1964

Sheet 1 of 3

4,466,427

IP Landscape Overview:

- Thomas Traction → Public Domain (centuries old)
- Improvements → patented
- Expired & fragmented IP across industries → High FTO

Predicate Analysis

Innovation

Tension Application: Free weights → Winch Mechanism

Portability: Fixed bed frame → detached fixture


Force Assessment: # Weights → Inline force gauge

Predicate Technology

Thomas Skeletal Traction
(Splint + Steinmann pin)

Novel Technology

Portable Skeletal Traction





Regulatory Pathway

FDA Regulatory Classification

Class 2 → 510(k) clearance pathway

External Traction Fixture with Invasive Pin

Substantial Equivalence → Thomas Traction

Relevant Regulatory Standards

ISO 13485 & ISO 14971 – System Quality/Risk

ISO 10993 – Biocompatibility

ASTM F1541 – Skeletal fixation testing

Clinical Testing

Future product differentiation/marketing

De novo pathway → not a target (cost)

Test against Thomas Traction → Equivalence

Product Codes

Thomas Traction → Grandfathered Device

HSP – Noninvasive traction splints (Fixture)

KTT – Invasive Skeletal fixation (Bone pin)

Business Model Canvas

Execution

R&D + Mechanical testing
Small Scale Clinical Trials
FDA 510(k) clearance

Revenue

1st: Device sale
2nd: Disposable kits + service
Optional: Leasing

Customers

Buyers: Hospitals
Users: Surgeons, staff
Stakeholders: Nurses, radiology

Costs

Upfront: R&D, regulatory
Ongoing: manufacturing, sales

Go-To-Market

Direct sales + distributors
Training + conferences
Clinical validation

Value

Portable traction
Constant force during movement
Imaging compatibility



Commercialization Strategy



Supply Chain

- Pin designed (in house) → device compatible
- Outsource component manufacturing
- Dual source → supply chain stability

Clinical Benefit Validating Cost

- No weights → reduce failure mode risk
- Imaging compatibility → better surgical outcomes
- Inline force assessment → responsive care
- Portable design → easy mobility across hospital
- Reusable design → cheaper over time for hospital

Unit Cost (COGS)	Device	\$803
	Traction Pin	\$20
Unit List Price	Device	\$1800
	Traction Pin	\$80
Gross Margin	Device	55.4%
	Traction Pins	75%

Summary

- Single sale (device) → Reoccurring sale (pins)
- Strong margins → support R&D + growth
- Competitive pricing with multifaceted clinical value
- Positioned for clean market entry and adaptation



Next Steps

Winch Design

Update winch design to add mechanism for limiting maximum allowable force and addressing human factors

Humans Factors Considerations

Improve crank arm-worm shaft interface, and provide large vs small force adjustment options.

Traction Bow and Pin Connection

Consider more details for design of traction bow and pin connection

Testing

Perform imaging and sterility verification, and conduct validation in clinical setting

Additional thanks to:

Our professors

Dr. Mario Cabodi
Professor Megan Lee
Professor Marc Gillette

Our industry mentors

Sam Grossman
Matt Sucec

Questions?

